

# PRDM11 Rabbit pAb

Catalog No.: A8502

## Basic Information

### Observed MW

58kDa

### Calculated MW

53kDa/57kDa

### Category

Primary antibody

### Applications

WB

### Cross-Reactivity

Mouse

## Background

May be involved in transcription regulation.

## Recommended Dilutions

WB 1:500 - 1:2000

## Immunogen Information

### Gene ID

56981

### Swiss Prot

Q9NQV5

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 228-477 of human PRDM11 (NP\_001243625.1).

### Synonyms

PRDM11;PFM8

## Contact

 | [www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification


Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

# Validation Data



Western blot analysis of extracts of various cell lines, using PRDM11 antibody (A8502) at 1:1000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.  
Lysates/proteins: 25ug per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Enhanced Kit (RM00021).  
Exposure time: 90s.