

# KIAA0513 Rabbit pAb

Catalog No.: A14647

## Basic Information

---

### Observed MW

47kDa

### Calculated MW

33kDa/45kDa/46kDa

### Category

Primary antibody

### Applications

WB

### Cross-Reactivity

Human, Mouse

## Background

---

## Recommended Dilutions

---

**WB** 1:500 - 1:2000

## Immunogen Information

---

### Gene ID

9764

### Swiss Prot

O60268

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-301 of human KIAA0513 (NP\_001273494.1).

### Synonyms

KIAA0513

## Contact

---

 | [www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification


### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thiomersal, 50% glycerol, pH7.3.

## Validation Data

---



Western blot analysis of extracts of various cell lines, using KIAA0513 antibody (A14647) at 1:1000 dilution.

Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25ug per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 1s.