

ZFP64 Rabbit pAb

Catalog No.: A15863

Basic Information

Observed MW

75kDa

Calculated MW

46kDa/48kDa/68kDa/72kDa/74kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Human, Mouse

Background

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

55734

Swiss Prot

Q9NPA5,Q9NTW7

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-170 of human ZFP64 (NP_060667.2).

Synonyms

ZFP64;ZNF338

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thiomersal, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using ZFP64 antibody (A15863) at 1:1000 dilution.

Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25ug per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 10s.