

TIMM50 Rabbit pAb

Catalog No.: A9991

Basic Information

Observed MW

39kDa

Calculated MW

27kDa/39kDa/50kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Human, Mouse

Background

Recommended Dilutions

WB 1:1000 - 1:4000

Immunogen Information

Gene ID

92609

Swiss Prot

Q3ZCQ8

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 190-456 of human TIMM50 (NP_001001563.1).

Synonyms

TIMM50;TIM50;TIM50L

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using TIMM50 antibody (A9991) at 1:1000 dilution.

Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25ug per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 1s.