

Phospho-Cytokeratin 18 (KRT18)-S33 Rabbit pAb

Catalog No.: AP0263 **1 Publications**

Basic Information

Observed MW

48kDa

Calculated MW

48kDa

Category

Primary antibody

Applications

WB,IHC

Cross-Reactivity

Human

Background

KRT18 encodes the type I intermediate filament chain keratin 18. Keratin 18, together with its filament partner keratin 8, are perhaps the most commonly found members of the intermediate filament gene family. They are expressed in single layer epithelial tissues of the body. Mutations in this gene have been linked to cryptogenic cirrhosis. Two transcript variants encoding the same protein have been found for this gene.

Recommended Dilutions

WB	1:500 - 1:1000
IHC	1:50 - 1:100

Immunogen Information

Gene ID	Swiss Prot
3875	P05783

Immunogen

A phospho specific peptide corresponding to residues surrounding S33 of human KRT18

Synonyms

KRT18;CK-18;CYK18;K18;keratin 18

Contact

 | www.abclonal.com


Product Information

Source	Isotype	Purification
Rabbit	IgG	Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Validation Data



Western blot analysis of extracts from HT29 cells untreated or treated with EGF using Phospho-Cytokeratin 18 (KRT18)-S33 antibody (AP0263).

Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25ug per lane.

Blocking buffer: 3% BSA.