

EIF3K Rabbit pAb

Catalog No.: A9969 **1 Publications**

Basic Information

Observed MW

25kDa

Calculated MW

24kDa/25kDa

Category

Primary antibody

Applications

WB, IHC

Cross-Reactivity

Human, Mouse, Rat

Background

Recommended Dilutions

WB	1:500 - 1:2000
IHC	1:100 - 1:200

Immunogen Information

Gene ID	Swiss Prot
27335	Q9UBQ5

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-218 of human EIF3K (NP_037366.1).

Synonyms

EIF3K;ARG134;EIF3-p28;EIF3S12;HSPC029;M9;MSTP001;PLAC-24;PLAC24;PRO1474;PTD001

Contact

 | www.abclonal.com


Product Information

Source	Isotype	Purification
Rabbit	IgG	Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using EIF3K antibody (A9969) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25ug per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 90s.