

# FAM234A Rabbit pAb

Catalog No.: A16580

## Basic Information

---

**Observed MW**  
67kDa

**Calculated MW**  
58kDa/59kDa

**Category**  
Primary antibody

**Applications**  
WB

**Cross-Reactivity**  
Human

## Background

---

Located in cell surface.

## Recommended Dilutions

---

**WB** 1:500 - 1:2000

## Immunogen Information

---

Gene ID	Swiss Prot
83986	Q9H0X4

**Immunogen**  
Recombinant fusion protein containing a sequence corresponding to amino acids 70-400 of human FAM234A (NP\_114428.1).

**Synonyms**  
FAM234A;C16orf9;ITFG3;gs19

## Contact

---

 | [www.abclonal.com](http://www.abclonal.com)

## Product Information


---

Source	Isotype	Purification
Rabbit	IgG	Affinity purification

**Storage**  
Store at -20°C. Avoid freeze / thaw cycles.  
Buffer: PBS with 0.01% thiomersal, 50% glycerol, pH7.3.

## Validation Data

---



Western blot analysis of extracts of various cell lines, using FAM234A antibody (A16580) at 1:1000 dilution.

Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25ug per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 3s.